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

Ver  mm

**Sample stage**

Ver	ThetaY
By use of channel <input type="text" value="CH4"/>	By use of channel <input type="text" value="CH4"/>
Mode <input type="radio"/> ABS <input checked="" type="radio"/> REL	Mode <input type="radio"/> ABS <input checked="" type="radio"/> REL
Start (mm) <input type="text" value="1"/>	Start (mm) <input type="text" value="0.8"/>
End (mm) <input type="text" value="-1"/>	End (mm) <input type="text" value="-0.8"/>
Step (mm) <input type="text" value="0.05"/>	Step (mm) <input type="text" value="0.02"/>
Integ (sec) <input type="text" value="0.1"/>	Integ (sec) <input type="text" value="0.1"/>
<b>Final position</b>	<b>Final position</b>
1st move to <input type="text" value="Differential peak"/>	<input checked="" type="radio"/> Move to Peak
2nd move to <input type="text" value="Original position"/>	<input type="radio"/> Move to Gravity

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The image shows a software window titled "GI stage alignment configuration". The window has a light blue background and a white border. At the top right, there are standard window control buttons (minimize, maximize, close). The main content area is divided into several sections:

- PD**: A section with a "Reset" button on the right. Below it is the text "Beam Stopper (Absolute)" and a text input field labeled "Ver" containing the value "9.482 mm".
- Sample stage**: A section containing two columns of controls.
  - Left Column:**
    - Ver**: "By use of channel" dropdown set to "CH4".
    - Mode**: Radio buttons for "ABS" and "REL", with "REL" selected.
    - Start (mm)**: Text input field with value "1".
    - End (mm)**: Text input field with value "-1".
    - Step (mm)**: Text input field with value "0.05".
    - Integ (sec)**: Text input field with value "0.1".
    - Final position**: Two dropdown menus. "1st move to" is set to "Differential peak" and "2nd move to" is set to "Original position".
  - Right Column:**
    - ThetaY**: "By use of channel" dropdown set to "CH4".
    - Mode**: Radio buttons for "ABS" and "REL", with "REL" selected.
    - Start (mm)**: Text input field with value "0.8".
    - End (mm)**: Text input field with value "-0.8".
    - Step (mm)**: Text input field with value "0.02".
    - Integ (sec)**: Text input field with value "0.1".
    - Final position**: Radio buttons for "Move to Peak" (selected) and "Move to Gravity".

At the bottom of the dialog, there are two buttons: "Update" and "Close". The "Update" button is highlighted with a yellow rectangular border. Below the buttons, a status bar contains the text "Input Ok. No changed value." and a small icon of three dots on the right.

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